

Name _____

Micro-slide Viewer Practice




Number of set _____

Using a Micro-Slide Viewer: Face the viewer so that as much light as possible falls on the white stage. Insert the numbered end of the slide holder into the slide slot of your viewer, moving it from your right to left. View with your eye close to the eye piece. With slide no. 1 in place, focus by turning the focus knob. Draw a **detailed** picture of what you see. Label it.

No. 1 Power _____

No. 2 Power _____



No. 3 Power _____



No. 4 Power _____

No. 5 Power _____

No. 6 Power _____



No. 7 Power _____

No. 8 Power _____

Extra Cr Power _____

